













Menus and keyboard Accelerators

eva

Menu	Keyboard Accelerator	Toolbar	Action
File New	Ctrl N		Create a new document
File Open	Ctrl O		Open an existing document
File Save	Ctrl S		Save the active document
File Information	Ctrl I		Patient information
File Start Recorder	R		Start real time recorder
File Playback	P		Start or stop playback
File Print	Ctrl P		Print the active document
Edit Copy	Ctrl C		Copy the selected graph or array and put it in the clipboard
View Zoom out	-		More general View (horizontal axis)
View Zoom in	+		More precise View (horizontal axis)
View Full Page	*		Full page view (horizontal axis)
View Move Left	←		Move the view to the left
View Move Right	→		Move the view to the right
View Scale out	Ctrl -		More general View (vertical axis)
View Scale in	Ctrl +		More precise View (vertical axis)
View Full Scale	Ctrl *		Automatic vertical scale
View Move Up	↑		Move up the selected graph
View Move down	↓		Move down the selected graph
Action Measures	M		Call or hide the measure cursor
Tools Options	O		Software options